

Citation for published version:

Mcclymont, D 2011, 'Multi-Step TPSR/QMS Technique to Study the Kinetics of NH₃-SCR', International Symposium on Nitrogen Oxides' Emission Abatement 2011, 4/09/11.

Publication date:

2011

Document Version

Early version, also known as pre-print

[Link to publication](#)

University of Bath

Alternative formats

If you require this document in an alternative format, please contact:
openaccess@bath.ac.uk

General rights

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

Take down policy

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

Multi-Step TPSR/QMS Technique to Study the Kinetics of NH₃-SCR

David W. J. McClymont,^a Joanna Ochońska,^{b,c} Stan T. Kolaczowski,^a Joanna Łojewska^c

^aDoctoral Training Centre, Centre for Sustainable Chemical Technologies, University of Bath, BA2 7AY, UK.

^bInstitute of Chemical Engineering of the Polish Academy of Sciences, Bałtycka 5, 44-100 Gliwice, Poland.

^cJagiellonian University, Faculty of Chemistry, Ingardena 3, 30-060 Kraków, Poland.

E-mail: D.W.J.McClymont@bath.ac.uk URL: <http://www.bath.ac.uk/csct>

1. What is NO_x?

- Nitric oxides are reactive gases, primarily NO (>90 %) and NO₂.
- They are pollutants involved in many atmospheric processes e.g. formation of photochemical smog and acid rain.
- They are produced as a result of high temperatures during the combustion of fuels.
- Legislation is in place to reduce NO_x emissions i.e. the European Waste Incineration Directive (WID) regulates activities that involve burning or gasification of waste(**Figure 1**).




Figure 1. Schematic of proposed biogas engine exhaust treatment system

4. TPSR/TPD Studies

- Coupled with Temperature Programmed Desorption (TPD), Temperature-Programmed Surface Reactions (TPSR) can be used to gain new insights into the mechanism of NO reduction on copper-exchanged zeolites. Results may also shed light on the active centres of the reaction and relative reactivities of surface and gaseous species.
- TPD is used to study the adsorption/desorption of NH₃ and NO (**Figure 4**).
- TPSR can be used to investigate reactions between adsorbed and gaseous species (**Figure 5**).




Figure 4. NO-TPD and NH₃-TPD plots of 3 copper exchanged zeolite powders.

REFERENCES



¹Kołodziej, A. & Łojewska, J. 2009. *Catalysis Today*, 147, S120-S124.

²Ochońska, J. *et al.* 2011 "Copper exchanged zeolite in structured reactor based on wire gauzes for ammonia reduction of nitrogen oxides from biogas turbines", NOEA 2011.

³Jodłowski, P. *et al.* 2011 "Modeling of structured reactor based on wire gauzes and zeolite catalyst for ammonia reduction of NOx biogas turbine", NOEA 2011

2. DeNO_x Process

- NH₃-Selective Catalytic Reduction (SCR) is an efficient, established method for NO_x removal. The desired reactions are:



- BUT there are some disadvantages including:

- Ammonia slip
- Size of the installation
- Thermal deactivation




Figure 2. Metal sheet catalyst supports

- Structured reactors based on metallic short channel structures (**Figure 2**) demonstrate improved mass and heat transfer properties [1] and can remedy these issues [2].

3. Catlab

- Hidden Analytical CATLAB-PCS microreactor (**Figure 3**) with integrated quadrupole mass spectrometer (QMS) system.
- QMS is a convenient and powerful method for analysing complex reaction mixtures. However, careful calibration is required to obtain reliable, quantitative data.
- There are a number of factors to consider when optimizing the analysis system for a particular application, including: detector selection, ion fragment patterns and relative sensitivity values.



Figure 3. Hidden Analytical CATLAB-PCS microreactor




Figure 5. Example NO-TPSR plot obtained from Cu-Y zeolite powder.

6. Conclusions

- Initial results suggest a Eley-Rideal-type mechanism i.e. adsorbed NH₃ reacts with gaseous NO on Cu⁺ active centres in the zeolite cavities.
- Obtained data was used for kinetic modelling [3].